EN Hydrogen Peroxide Solution with neutralizing tablets

For cleaning and disinfecting, neutralizing, storage and inserting.
For all types of contact lenses.

Composition hydrogen peroxide solution:

HYDROGEN PEROXIDE, PHENACETINUM CRYSTALLINUM, SODIUM CHLORIDE, PVP K30 EP, PURIFIED WATER.

Composition neutralizing tablet:

BORIC ACID, POTASSIUM PENTABORATE, POLYVINYL ALCOHOL, VITAMINE B12, MICROBIAL CATALASE.

Precautions and instructions for use:

The hydrogen peroxide solution should always be used in combination with the enclosed neutralizing tablets and a lens case with air valve.
Wash your hands carefully before removing the lenses from your eyes or the lens case. Dry your hands using a clean towel.
Do not use when tamper-evident seal is broken before the first time use. Unscrewing the cap from the bottle will break the tamper-evident seal.
Rinse lens case (with air valve) with saline solution before filling with hydrogen peroxide solution, to remove all remnant solution.
Place the left contact lens in the cup marked 'L' and close the cup. Place the right contact lens in the cup marked 'R' and close the cup.
Fill the lens case with hydrogen peroxide solution up to the fill line.
Add one neutralizing tablet to the lens case filled with hydrogen peroxide.
The solution will turn pink.
Allow the lenses to neutralize for a minimum of 6 hours. After neutralizing the lenses may be inserted.
After treatment the lens can appear slightly pinkish, this color will soon disappear during wearing.
Repeat the process after every 24 hours of storage.
Do not touch the dropper tip of the bottle as this may contaminate the solution.
Close the bottle tightly immediately after use.
Replace your lens case with air valve frequently to prevent infection due to contamination.

Storage conditions:

Use within 2 months after opening.
The tablets are for single use only.
Store between 4 °C and 25 °C.
Protect from sunlight.
Keep out of reach of children (see image 1).
Keep the solution in its original packaging.

Warnings and contra-indications:

Do not insert your contact lenses directly in your eye without neutralizing your contact lenses since this can lead to pain and burning eyes.
Do not use if you are allergic to one of the ingredients.
Always use a lens case with air valve.
Do not insert the lenses in the case the solution has not turned pink, repeat the neutralizing process.
Stop wearing contact lenses when you experience eye irritation or when you have an eye infection and consult a doctor.
If you are prescribed eye medication, consult your ophthalmologist or contact lens specialist before using the medication in combination with contact lens solution.
Do not use when the packaging is damaged (see image 2).
Do not use after the expiry date.
Do not mix your contact lens care system without consulting your contact lens specialist.
Do not change or use with other contact lens solution, unless it is indicated in the instructions for use.
Never use tap water as an alternative for contact lens solution.
Stop using the product and inform the legal manufacturer and the national health authorities when you experience a serious health issue in relation with this medical device.

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