



# MONOSEPT



3% hydrogen peroxide system with neutralisation disc for soft contact lenses, including silicone hydrogel lenses.

## PURPOSE

Disinfecting, removing protein deposits and storing all soft contact lenses.

## BENEFITS

- Preservative-free
- Thorough disinfection thanks to a 3% hydrogen peroxide solution
- Disinfection cycle takes only 6 hours

## COMPOSITION

3% hydrogen peroxide solution, Phosphonic acid stabiliser, NaCl, Phosphate buffers, purified water

## AVAILABLE IN

- **120 ml Pack** (1 x 120 ml + lens case with disc)  
EAN - 8717185944761
- **360 ml Pack** (1 x 360 ml + lens case with disc)  
EAN - 8717185944778
- **Multi Pack** (5 x 360 ml + 5 lens cases with disc)  
EAN - 8717185940343